Principes fondamentaux de la tDCS

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
1. Koenigs, M., & Grafman, J. (2009). The functional neuroanatomy of depression: Distinct roles for ventromedial and dorsolateral prefrontal cortex. *Behavioural Brain Research*, 201(2), 239. <https://doi.org/10.1016/j.bbr.2009.03.004>
2. Campbell S et al. (2004). Lower hippocampal volume in patients suffering from depression: a meta-analysis. *Am J Psychiatry*, 161(4):598-607. <https://doi.org/10.1176/appi.ajp.161.4.598>
3. Hamilton, J. P., Siemer, M., & Gotlib, I. H. (2008). Amygdala volume in Major Depressive Disorder: A meta-analysis of magnetic resonance imaging studies. *Molecular Psychiatry*, 13(11), 993. <https://doi.org/10.1038/mp.2008.57>
4. Lefaucheur JP et al. (2016). Evidence-based guidelines on the therapeutic use of transcranial direct current stimulation (tDCS). *Clinical Neurophysiology*, 128(1), 56-92. <https://doi.org/10.1016/j.clinph.2016.10.087>
5. Nitsche, M. A., & Paulus, W. (2000). Excitability changes induced in the human motor cortex by weak transcranial direct current stimulation. *The Journal of Physiology*, 527(Pt 3), 633-639. <https://doi.org/10.1111/j.1469-7793.2000.t01-1-00633.x>

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1. Clayden, J. D. (2013). Imaging connectivity: MRI and the structural networks of the brain. *Functional Neurology*, 28(3), 197. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3812744/>


2. Adams, T. G., Cisler, J. M., Kelmendi, B., George, J. R., Kichuk, S. A., Averill, C. L., Anticevic, A., Abdallah, C. G., & Pittenger, C. (2021). Transcranial direct current stimulation targeting the medial prefrontal cortex modulates functional connectivity and enhances safety learning in obsessive-compulsive disorder: Results from two pilot studies. *Depression and Anxiety*, 39(1), 37-48. <https://doi.org/10.1002/da.23212>

3. Antonenko, D., Fromm, A. E., Thams, F., Grittner, U., Meinzer, M., & Flöel, A. (2023). Microstructural and functional plasticity following repeated brain stimulation during cognitive training in older adults. *Nature Communications*, 14(1), 1-13. <https://doi.org/10.1038/s41467-023-38910-x>

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
1. Fregni et al. (2021). Evidence-Based Guidelines and Secondary Meta-Analysis for the Use of Transcranial Direct Current Stimulation in Neurological and Psychiatric Disorders. International Journal of Neuropsychopharmacology, 24(4), 256–313. <https://doi.org/10.1093/ijnp/pyaa051>
2. Nikolin et al. (2023). Time-course of the tDCS antidepressant effect: An individual participant data meta-analysis. Progress in Neuro-Psychopharmacology and Biological Psychiatry, 125, 110752. <https://doi.org/10.1016/j.pnpbp.2023.110752>
3. Løokene et al. (2022). Reduction of symptoms in patients with major depressive disorder after transcranial direct current stimulation treatment: A real-world study. Journal of Affective Disorders Reports, 8, 100347. <https://doi.org/10.1016/j.jadr.2022.100347>
4. Bikson et al. (2016). Safety of Transcranial Direct Current Stimulation: Evidence Based Update 2016. Brain Stimulation, 9(5), 641–661. <https://doi.org/10.1016/j.brs.2016.06.004>

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Woodham, R.D., Selvaraj, S., Lajmi, N., et al. (2025).

Home-based transcranial direct current stimulation treatment for major depressive disorder: a fully remote phase 2 randomized sham-controlled trial.

Nature Medicine, 31, 87–95. <https://doi.org/10.1038/s41591-024-03305-y>

Oh, J., Jang, I., Jeon, S., & Chae, H. (2022).

Effect of Self-administered Transcranial Direct Stimulation in Patients with Major Depressive Disorder: A Randomized, Single blinded Clinical Trial.

Clinical Psychopharmacology and Neuroscience, 20(1), 87–96. <https://doi.org/10.9758/cpn.2022.20.1.87>

Zhou Q, Yu C, Yu H, Zhang Y, Liu Z, Hu Z, Yuan TF, Zhou D. (2020).

The effects of repeated transcranial direct current stimulation on sleep quality and depression symptoms in patients with major depression and insomnia.

Sleep Med. Jun;70:17–26. doi: 10.1016/j.sleep.2020.02.003.

Sharafi E, Taghva A, Arbabi M, Dadarkhah A, Ghaderi J. (2019).

Transcranial Direct Current Stimulation for Treatment-Resistant Major Depression: A Double-Blind Randomized Sham Controlled Trial.

Clin EEG Neurosci. Nov;50(6):375–382. doi: 10.1177/1550059419863209.

Salehinejad, M. A., Ghanavai, E., Rostami, R., & Nejati, V. (2017).

Cognitive control dysfunction in emotion dysregulation and psychopathology of major depression (MD): Evidence from transcranial brain stimulation of the dorsolateral prefrontal cortex (DLPFC).

Journal of affective disorders, 210, 241–248. <https://doi.org/10.1016/j.jad.2016.12.036>

Études notables

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Valiengo, L., Benseñor, I. M., Goulart, A. C., Zanao, T. A., Boggio, P. S., Lotufo, P. A., Fregni, F., & Brunoni, A. R. (2013). **The Sertraline Versus Electrical Current Therapy For Treating Depression Clinical Study (Select-Tdcs): Results Of The Crossover And Follow-Up Phases.** Depression and Anxiety, 30(7), 646–653. <https://doi.org/10.1002/da.22079>

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Valiengo LC, Goulart AC, de Oliveira JF, Benseñor IM, Lotufo PA, Brunoni AR. (2017). **Transcranial direct current stimulation for the treatment of post-stroke depression: results from a randomised, sham-controlled, double-blinded trial.** J Neurol Neurosurg Psychiatry. Feb;88(2):170–175. doi: 10.1136/jnnp-2016-314075.

tDCS vs Escitalopram

Brunoni, A. R., Moffa, A. H., Sampaio-Junior, B., Borriane, L., Moreno, M. L., Fernandes, R. A., Veronezi, B. P., Nogueira, B. S., Aparicio, L. V. M., Razza, L. B., Chamorro, R., Tort, L. C., Fraguas, R., Lotufo, P. A., Gattaz, W. F., Fregni, F., Benseñor, I. M., & ELECT-TDCS Investigators (2017). **Trial of Electrical Direct-Current Therapy versus Escitalopram for Depression.** The New England journal of medicine, 376(26), 2523–2533. <https://doi.org/10.1056/NEJMod1612999>

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Vigod SN, Murphy KE, Dennis CL, Oberlander TF, Ray JG, Daskalakis ZJ, Blumberger DM. (2019). **Transcranial direct current stimulation (tDCS) for depression in pregnancy: A pilot randomized controlled trial.** Brain Stimul. Nov-Dec;12(6):1475–1483. doi: 10.1016/j.brs.2019.06.019.